



Email : sales@laboratoryinstrumentindia.com

Product Name :

Gel, ultrasound, bottle, min. 250ml

Product Code :

LBNY-0021-22000126



www.laboratoryinstrumentindia.com

Description :

A non-sterile medium designed to be applied to a patient's unbroken skin surface to provide a coupling between an analytical device (e.g., ultrasound transducer, optical glucose monitoring system) and the patient, allowing for the emission and reception of signals (e.g., ultrasound, light) that pass through the skin during an examination. It may be used by a healthcare professional or layperson. After application, this device cannot be reused. It is in the form of a non-sticky, fluid-like gel that may also assist the operator's ability to move the analytical device smoothly over the skin to better analyse the site of interest.

Technical Specification :

Excellent transmission of ultrasound waves for a broad range of frequencies.
Reduces or eliminates static energy.
Excellent gliding properties.
Sufficient viscosity to avoid dripping.
Hypoallergenic.
Bottle size between 250 and 260 ml.
Does not stain fabrics or damages equipment.
Non-sterile product.

Laboratory Instrument India

