



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Needle,scalp vein,25G,ster,disp

**Product Code :**

LBNY-0021-22000128



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Material: Sharp needle, stainless steel, siliconised, for atraumatic insertion.

Disposable.

Sterile.

External diameter expressed in Gauge and mm.

Colour code/external diameter: Visible on the wings.

Components: Needle with protecting cap, wings, tube including Luer lock connector with cap (all parts fit together to form a unit).

Length expressed in mm.

Initial sterilisation method: Ethylene oxide gas.

Soft flexible tubing, PUR (polyurethane) or PVC, length approx. 150-300mm, including Luer connector with cap.

Butterfly wings, PVC, for better handling and fixation.

Size selected: 25G (0.50 x 20mm) orange.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.043cdm

Estimated weight: 0.006kg

---

# Laboratory Instrument India

