



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Needle,scalp vein,21G,ster,disp

**Product Code :**

LBNY-0021-22000129



www.laboratoryinstrumentindia.com

**Description :**

External diameter expressed in Gauge and mm

Material: Sharp needle, stainless steel, siliconised, for atraumatic insertion.

Butterfly wings, PVC, for better handling and fixation.

Disposable.

Sterile.

Initial sterilisation method: Ethylene oxide gas.

Soft flexible tubing, PUR (polyurethane) or PVC, length approx. 150-300mm, including Luer connector with cap.

Size selected: 21G (0.80 x 20mm) green.

Length expressed in mm

Colour code/external diameter: Visible on the wings.

Components: Needle with protecting cap, wings, tube including Luer lock connector with cap (all parts fit together to form a unit).

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.038cdm

Estimated weight: 0.005kg

---

# Laboratory Instrument India

