



Email : sales@laboratoryinstrumentindia.com

Product Name :
Argon Gas Cylinder

Product Code :
LBNY-0004-820046



Description :

Argon Gas Cylinder

Technical Specification :

Argon gas cylinder
Argon gas
40l Gas Tank

Pure argon is colourless, tasteless, odourless, inert to all materials at all temperature and pressure. Argon is heavier than air and hence will collect in low lying areas such as ducts and drains. Its an asphyxiant. Argon is available in multiple packages in purity level of 99.9% or more. Argon is used:

Colourless, tasteless, odourless and inert to all materials at all temperature and pressure

Shielding gas for MIG and TIG welding processes

Welding and cutting a wide range of ferrous and non-ferrous materials

Plasma cutting



Laboratory Instrument India

