



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Argon Gas Cylinder

**Product Code :**  
LBNY-0004-820046



**Description :**

Argon Gas Cylinder

**Technical Specification :**

Argon gas cylinder  
Argon gas  
40l Gas Tank

Pure argon is colourless, tasteless, odourless, inert to all materials at all temperature and pressure. Argon is heavier than air and hence will collect in low lying areas such as ducts and drains. Its an asphyxiant. Argon is available in multiple packages in purity level of 99.9% or more. Argon is used:

Colourless, tasteless, odourless and inert to all materials at all temperature and pressure

Shielding gas for MIG and TIG welding processes

Welding and cutting a wide range of ferrous and non-ferrous materials

Plasma cutting



**Laboratory Instrument India**

