



Email : sales@laboratoryinstrumentindia.com

Product Name :
Gas Cylinder Trolley

Product Code :
LBNY-0004-820047



Description :

Gas Cylinder Trolley

Technical Specification :

Gas cylinder trolley

Made from best quality of mild steel material.

These Cylinder Truck are effective to carry different types of Gas Cylinder Trolley from one place to other.

These Cylinder Truck have 1 swivel castors wheels with double ball bearing that assist easy transportation of materials in required directions.

These are widely used in oxygen bottle mfg. plant, factories, engineering company, assemble department, etc.

Gas cylinders: 50 litres(each: 2)

Ø gas cylinder, max, (mm): 250

Material: Steel

Surface: Galvanized

Tyres: Solid rubber

Total payload (kg): 200



Laboratory Instrument India

