



Email : sales@laboratoryinstrumentindia.com

Product Name :

Uterine Balloon Tamponade, Sterile, S.U.

Product Code :

LBNY-0021-22000154



www.laboratoryinstrumentindia.com

Description :

The Uterine Balloon Tamponade (UBT) is a sterile, single-use medical device designed to control postpartum hemorrhage (PPH) by applying pressure to the uterus, thereby reducing blood loss and stabilizing the woman undergoing the condition. The balloon is inserted into the uterus, filled with saline solution or sterile water, and effective tamponade occurs rapidly to stop the bleeding.

Uterine balloon tamponades (UBT) are medical devices that apply pressure to the uterus to reduce blood loss, stabilize the woman undergoing postpartum haemorrhage (PPH), and reduce the need for surgical interventions and blood transfusions.

Technical Specification :

Material: medical grade elastomer

Capacity: 500 ml – 700 ml filling capacity before balloon stretch

Stretch capacity: minimum 2000 ml

Shape: Diameter 80 mm and cylindrical shape

Thickness: Thin walled for low expansion pressure requirement



Laboratory Instrument India

