



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Sachet,tablet,plastic,10x16cm/PAC-100

**Product Code :**

LBNY-0021-22000155



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Sachet, transparent, virgin polyethylene, self sealing with interlock type closure on 1 short side, securely sealed along the 2 long sides.

Plastic Tablet bag, sealable with interlock type closure, 10 x 16cm, pack of 100.

**Technical Specification :**

Presented with three white write-on horizontal stripes, clearly and cleanly printed on one side of the sachet, of approximate dimensions: 100mm across x 15mm in width. Dimensions: active length (length below the base of the seal) of 150mm, length above the top of the seal of 10mm, width of 100mm and thickness of 0.05mm (50 microns).



**Laboratory Instrument India**