



Email : sales@laboratoryinstrumentindia.com

Product Name :

Sachet,tablet,plastic,10x16cm/PAC-100

Product Code :

LBNY-0021-22000155



www.laboratoryinstrumentindia.com

Description :

Sachet, transparent, virgin polyethylene, self sealing with interlock type closure on 1 short side, securely sealed along the 2 long sides.

Plastic Tablet bag, sealable with interlock type closure, 10 x 16cm, pack of 100.

Technical Specification :

Presented with three white write-on horizontal stripes, clearly and cleanly printed on one side of the sachet, of approximate dimensions: 100mm across x 15mm in width. Dimensions: active length (length below the base of the seal) of 150mm, length above the top of the seal of 10mm, width of 100mm and thickness of 0.05mm (50 microns).



Laboratory Instrument India