



Email : sales@laboratoryinstrumentindia.com

Product Name :
Coverall, protection, CatIII, type 6b, XL

Product Code :
LBNY-0021-22000158



www.laboratoryinstrumentindia.com

Description :

Spray/aerosol-penetration resistant, biohazard-protective coverall, for use in isolation units for infection prevention and control against viral penetration.

Technical Specification :

Elasticated hood around face.
Elasticated cuffs and ankles.
Material: Lightweight, do not contain rubber/ latex.
Antistatic treated on both sides.
Sleeves with elasticated thumb loop.
Stitched-in neck label indicating the type and performance of the suit against the below mentioned standards.
Non-sterile
Single Use, disposable
Color: White
Fabric is Infective agent tested against viral penetration at minimum 1.75kPa (class 2) (or equivalent international standard)
Stitched or Bound seams.
Zipper with re-sealable flap protecting leakage through seams.
Size: XL.

Laboratory Instrument India

