



Email : sales@laboratoryinstrumentindia.com

Product Name :
Emergency food ration,500g

Product Code :
LBNY-0022-220002



www.laboratoryinstrumentindia.com

Description :

Fortified and cereal based dry compressed food, 9 bars per pack of 500g, 24 packs per carton. Each bars weighs minimum 55g.

Technical Specification :

Ingredients

No ingredients of animal origin

Allergen statement: May contain traces of milk.

Baked wheat, partially hydrogenated vegetable fat, sugar, vegetable soya protein, malt, vitamins and minerals.



Laboratory Instrument India