



Email : sales@laboratoryinstrumentindia.com

Product Name :
Therap.milk CAN 400g

Product Code :
LBNY-0022-220003



www.laboratoryinstrumentindia.com

Description :

The liquid therapeutic diet is intended for stabilisation phase of treatment of children diagnosed with SAM. It is intended to stabilize the child's metabolism and support rehydration. The specified quantity of powdered therapeutic milk is mixed with the specified quantity of water boiled and cooled down to not below 70°C to obtain a defined quantity of liquid therapeutic diet with an energy density of approximately 75 kcal/ 100ml. It comes as powder packed in a canister to prepare a liquid therapeutic diet with an energy density of approximately 75kcal/100 ml. Therapeutic milk diet with added vegetable fat, carbohydrates, vitamins and minerals.

Technical Specification :

Ingredients: maltodextrin (carbohydrate), vegetable oils, saccharose, skim milk powder or other dairy powder, emulsifier (eg lecithin) vitamin and mineral premix.
Milk based white or pale yellowish fine powder; free from impurities, coloured particles, caking or lumps.



Laboratory Instrument India