



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
F-100 Therap.milk CAN 400g

**Product Code :**  
LBNY-0022-220004



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

It comes as powder packed in a canister to prepare a liquid therapeutic diet with an energy density of approximately 100 kcal/100ml. Therapeutic milk diet with added vegetable fat, carbohydrates, vitamins and minerals.

**Technical Specification :**

The F100 liquid therapeutic diet is intended to be used for the second phase (nutritional rehabilitation phase) of treatment of children diagnosed with SAM to achieve rapid weight gain. It must be administered under medical supervision. The specified quantity of powdered therapeutic diet is mixed with the specified quantity of boiled and cooled down to not below 70°C water to obtain defined quantity of liquid therapeutic diet with an energy density of approximately 100 kcal/100ml.



**Laboratory Instrument India**