



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Therapeutic CMV,tin 800g/CAR-6

**Product Code :**  
LBNY-0022-220005



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Therapeutic complex of vitamins and minerals (CMV) is for a nutritional fortification for making therapeutic foods in the dietetic treatment of severe malnutrition. Use to prepared phase 1 milk (F-75) and enriched high energy milk (F-100). Therapeutic CMV is packed in 800g metal tins, in cartons of 6 tins. Each tin contains a measuring scoop to weigh 6.8 g.

**Technical Specification :**

It must only be used in the treatment of Severe Acute Malnutrition, therefore under medical supervision. The composition of therapeutic CMV complies with the latest recommendations for mineral and vitamin enrichment in the dietetic treatment of Severe Acute Malnutrition.



**Laboratory Instrument India**