



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
RUTF biscuit, packs, 510g

**Product Code :**  
LBNY-0022-220007



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

RUTF biscuit is to be used by humanitarian agencies, governmental and non-governmental organisations for the treatment of Severe Acute Malnutrition (SAM) in any cultural setting. RUTF biscuit may be used in a wide variation of climatic zones and may be the sole source of food, except for water and breast milk, during the period of use.

Ready-to-Use Therapeutic Food (RUTF) biscuit, a high-energy fortified food equivalent to that of the rehabilitation diet but with a higher energy density.

Each bar is 2 biscuits of 28.4 g.

Packed in cartons of 24 packs where each pack contains 9 bars.

**Technical Specification :**

Ready to use: be eaten directly (no cooking/mixing/dilution required) or crumbled into drinking water and eaten as porridge.

RUTF biscuits are compressed bars, manufactured from a mixture of cereal, milk powder, vegetable oil and carbohydrates, with added vitamins and minerals.

Storage conditions: as defined by the manufacturer, no refrigeration required.

Shelf life: 4 years

---

# Laboratory Instrument India

