



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Manual centrifuge

**Product Code :**  
LBNY-0022-220009



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Manual centrifuge to be used for preparing sample with portable fluorometer, vitamin A in foods and portable photometer, iron in foods.

**Technical Specification :**

Centrifuge can be clamped to a bench or a table  
Low maintenance required  
Self-lubricating gear.



**Laboratory Instrument India**