



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electrical Machine Trainer

Product Code :
LBNY-0002-620005



Description :

Electrical Machine Trainer

Technical Specification :

Facilitates easy and safe wiring by students due to use of 4mm sturdy banana patch cords and shrouded socket arrangements.

All machines are mounted on finely painted sturdy base frame with easy machine .

Use of gear coupling facilitates.

With due emphasis on student safety machines operate upto levels and upto 1500 RPM, without compromising on didactic use. Able to draw all graphs.

Electric Machine, DC, AC, Shunt, Squarrel Cage, Compound, Techo meter, Brake

Electrical Machine Trainer, Electrical Machine Trainer Equipments, Electrical Machine Trainer Tools, Electrical Machine Trainer Tool Kits, Electrical Machine Trainer Manufacturers, Electrical Machine Trainer Suppliers from India, China, Kenya



Laboratory Instrument India

