



Email : sales@laboratoryinstrumentindia.com

Product Name :
Reagents Vit A/Food

Product Code :
LBNY-0022-220011



www.laboratoryinstrumentindia.com

Description :

Kit with extraction vial, syringe and needle, used for the measurement with portable fluorometer to determine quantity of vitamin A content in foods, pack of 100.

Technical Specification :

Each extraction vial contains 2 ml of reagents, developed to extract and determine quantity of vitamin A from the sample.

Supplied as set containing:

100 x Reagent vial

Number of measurements: 100

100 x Hypodermic needle, 21G

100 x Single use syringe, 1ml

Composition of reagent:

Propan-2-ol

Ethanol

N-Hexane.



Laboratory Instrument India

