



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Supercereal Plus (CSB++)/BAG-1,5KG

**Product Code :**  
LBNY-0022-2200015



www.laboratoryinstrumentindia.com

**Description :**

SUPER CEREAL plus is prepared from heat treated maize and de-hulled soya beans, sugar, dried skim milk, refined soya bean oil, vitamins and minerals. Super Cereal Plus Corn Soya Blend/BAG-1,5kg. Fortified flour of corn and soy, vegetable oil, vitamins and minerals with added sugar and milk.

**Technical Specification :**

Ingredients / Formula

De-hulled soya beans: 20%  
Dried skim milk powder: 8%  
Refined soya bean oil: 3%  
Vitamin and Mineral premix: 0.20%  
Maize (white or yellow): 58.30%  
Sugar: 9%  
Potassium chloride: 0.27%  
Dicalcium phosphate anhydrous: 1.23%.



**Laboratory Instrument India**

