



Email : sales@laboratoryinstrumentindia.com

Product Name :
Supercereal Plus (CSB++)/BAG-1,5KG

Product Code :
LBNY-0022-2200015



www.laboratoryinstrumentindia.com

Description :

SUPER CEREAL plus is prepared from heat treated maize and de-hulled soya beans, sugar, dried skim milk, refined soya bean oil, vitamins and minerals. Super Cereal Plus Corn Soya Blend/BAG-1,5kg. Fortified flour of corn and soy, vegetable oil, vitamins and minerals with added sugar and milk.

Technical Specification :

Ingredients / Formula

De-hulled soya beans: 20%
Dried skim milk powder: 8%
Refined soya bean oil: 3%
Vitamin and Mineral premix: 0.20%
Maize (white or yellow): 58.30%
Sugar: 9%
Potassium chloride: 0.27%
Dicalcium phosphate anhydrous: 1.23%.



Laboratory Instrument India

