



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Supercereal Plus (WSB++)/BAG-1,5KG

**Product Code :**  
LBNY-0022-220016



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

SUPER CEREAL plus is prepared from heat treated wheat and de-hulled soya beans, sugar, dried skim milk, refined soya bean oil, vitamins and minerals. Super Cereal Plus Wheat Soya Blend/BAG-1,5kg. Fortified flour of wheat and soy, vegetable oil, vitamins and minerals with added sugar and milk.

**Technical Specification :**

Wheat 52.30%  
De-hulled soya beans: 25%  
Vitamin and Mineral premix: 0.20%  
Dicalcium phosphate anhydrous: 1.23%  
Dried skim milk powder: 8%  
Sugar: 9%  
Refined soya bean oil: 4%  
Potassium chloride: 0.27%.



**Laboratory Instrument India**

