



Email : sales@laboratoryinstrumentindia.com

Product Name :

Scale,beam type,child/adult,220kgx50g

Product Code :

LBNY-0022-2200020



www.laboratoryinstrumentindia.com

Description :

Beam type mechanical scale for adults

Display: easily readable in low light working situations, black coloured numbers on silver surface. Can be read from both sides.

With easy zero adjustment possibility.

With stabilizing mechanism for faster reading of results

Measuring range: from 5 up to a minimum of 180kg.

Readout in metric (kg)

Adjustable feet allow for horizontal levelling

Minimum graduation: 100g

Reading time in seconds: Hand pole attached.

No height measuring rod included.

Technical Specification :

Materials:

Body: metal only (stainless steel or epoxy coated steel).

Calibration screw: stainless steel

Body and mechanics are splash proof; shock resistant.

Design allows rough handling.

Mechanics and other vital parts are made of rust proof materials.

Ferrous parts shall be treated with a rust inhibitor against corrosion, plastic parts shall be UV resistant, preferably acrylonitrile butadiene styrene (ABS) or similar.

Smooth surface and finishing allow easy cleaning and disinfection.

Easy maintenance and repair in low tech settings possible.



Laboratory Instrument India