



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hand Pallet Truck

Product Code :
LBNY-0004-820051



Description :

Hand Pallet Truck

Technical Specification :

Hand Pallet Truck
Rubber coated handle
Epoxy Powder coated body
Hard chrome plated piston
Single piece pump with galvanized body
Large diameter polyurethane steering wheel
Tandem polyurethane load rollers
Foot Release / Hand Release
Entry slide for easy entry in closed pallets.
High load carrying capacity polyurethane wheels & rollers to protect floors.
Adjustable Push Rods.
4mm thick high tensile sheet for extra strenght of forks & body.
Push rods made from solid square bars instead of pipe.
Heavy duty additional angle supports on lifting lever frame.
Working load: 1,500- 1600 kg
Lift height for fork: (205-250) mm
Fork height: 85-100 mm
Fork length mm; 1,000-1200
Fork outer width mm: 550-600



Laboratory Instrument India