



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Scale,beam type,clinical,infant,16kgx5g

**Product Code :**

LBNY-0022-2200021



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Mechanical infant clinical scale, beam type, with tray, for use up to 16kg with 10g precision, with splash cover. Used for weighing infants in the context of growth monitoring.

Beam type mechanical scale for infants

Display: easily readable in low light working situations, signal white coloured numbers on black surface.

With easy zero adjustment possibility.

With stabilizing mechanism for faster reading of results.

Measuring range: 0-16kg

Minimum graduation: 5g

Readout in metric (kg)

Removable curved tray.

Reading time in seconds: about 30 seconds

Built in round level bubble with indication of centre.

Adjustable feet allow for horizontal levelling.

**Technical Specification :**

**Materials**

Body: metal only (stainless steel or epoxy coated steel).

Calibration screw: stainless steel

Body and mechanics must be splash proof, shock resistant and light-weight.

Design must allow rough handling.

Ferrous parts shall be treated with a rust inhibitor against corrosion ,

Plastic parts shall be UV resistant, preferably acrylonitrile butadiene styrene (ABS) or similar.

---

Smooth surface/finishing must allow easy cleaning/disinfection.  
Mechanics and other vital parts are made of rust proof materials  
Tray: powder-coated, light grey colour  
Easy maintenance and repair in low tech settings must be possible.



**Laboratory Instrument India**