



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
MUAC,Child 11.5 Red

**Product Code :**  
LBNY-0022-2200023



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

To measure the mid upper arm circumference for rapid screening of nutritional status of children. Target population: Malnourished children or those at risk, aged 6-59 months.

**Technical Specification :**

Tape for measuring child's Mid Upper Arm Circumference (MUAC).  
With three incisions to facilitate tape insertion and fixation, including one wider result window.  
Double-sided, featuring user instructions in English ('tape method') on the reverse side.  
Print: Permanent, resistant to solvents; easily readable in low light working situations.  
Measuring range up to 26.5 cm.  
Accuracy:  $\pm 1$  mm of the maximum measurement (26.5 cm).  
Working/storage temperature range: 0 - 50 °C.  
Graduated with 1 mm precision with thicker line at 21.0 cm.

**Material:**

Materials used are free from phthalates, heavy metals and other dangerous matters as per the latest Candidate List of substances of very high concern for Authorisation published on the website of the European Chemical Agency.  
Non-tear stretch-resistant polypropylene.  
Minimum thickness 0.3 mm.



**Laboratory Instrument India**