



Email : sales@laboratoryinstrumentindia.com

Product Name :
MUAC,Adult,w/o color code

Product Code :
LBNY-0022-220024



www.laboratoryinstrumentindia.com

Description :

Intended use:

To measure the mid upper arm circumference for rapid screening of nutritional status of pregnant/lactating women.

Adults' Mid Upper Arm Circumference (MUAC) measuring tape with a range of up to 50 cm. Pack of 50 tapes with pictorial and text instructions for use printed on the reverse side, reusable, in English language.

Technical Specification :

Tape for measuring adult's Mid Upper Arm Circumference (MUAC).

Double-sided, featuring user instructions in English ('tape method') on the reverse side.

Print: Permanent, resistant to solvents; easily readable in low light working situations.

With three incisions to facilitate tape insertion and fixation, including one wider result window.

Measuring range up to 50 cm.

Graduated with 1 mm precision.

Accuracy: ± 1 mm of the maximum measurement (50 cm).

Working/storage temperature range: 0 - 50 °C.

Material:

Non-tear stretch-resistant polypropylene.

Materials used are free from phthalates, heavy metals and other dangerous matters as per the latest Candidate List of substances of very high concern for Authorisation published on the website of the European Chemical Agency.

Minimum thickness 0.3mm.



Laboratory Instrument India