



Email : sales@laboratoryinstrumentindia.com

Product Name :
ID Bracelet,child,blue

Product Code :
LBNY-0022-2200025



www.laboratoryinstrumentindia.com

Description :

Closed with a strong, single use clip, that once clicked in to position cannot be removed.
Personal information to be written directly on the bracelet with common ballpoint pen, non-erasable, or in the paper strip then inserted inside the transparent bracelet.
Single use, tamper proof, transparent, non-stretch personal identification bracelet for paediatric use.
Colour: Blue.
Strong and durable, soft to the touch. Atraumatic and hypo-allergenic.
Dimensions: approximately 1 cm x 25 cm (w x l). Length: adjustable at closure, the remaining free end of the strap may be trimmed adjacent its emergence from the clip.

Technical Specification :

Material:
Identification paper (when applicable): writable with common ballpoint pen. Paper should be easy to insert into the accompanying plastic band.
Band: supple non-stretch plastic, LDPE-based, vinyl or equivalent. Resistant to water, common solvents and soaps. Single use.



Laboratory Instrument India

