



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
ID Bracelet,child,red

**Product Code :**  
LBNY-0022-2200026



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Strong and durable, soft to the touch. Atraumatic and hypo-allergenic.  
Closed with a strong, single use clip, that once clicked in to position cannot be removed.  
Personal information to be written directly on the bracelet with common ballpoint pen, non-erasable, or in the paper strip then inserted inside the transparent bracelet.  
Single use, tamper proof, transparent, non-stretch personal identification bracelet for paediatric use.  
Colour: red.  
Dimensions: approximately 1 cm x 25cm (w x l). Length: adjustable at closure, the remaining free end of the strap may be trimmed adjacent its emergence from the clip.

**Technical Specification :**

Material:  
Identification paper (when applicable): writable with common ballpoint pen. Paper should be easy to insert into the accompanying plastic band.  
Band: supple non-stretch plastic, LDPE-based, vinyl or equivalent. Resistant to water, common solvents and soaps. Single use.



**Laboratory Instrument India**

