



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Weighing trousers,infant,25kg

**Product Code :**  
LBNY-0022-2200027



www.laboratoryinstrumentindia.com

**Description :**

Weighing trousers for use with spring scales. One size, pack of 5 pieces in a carry bag.  
Used for weighing infants in the context of growth monitoring and for anthropometric surveys, for children aged 6-24 months.

**Technical Specification :**

Trousers:

Strong and easy maintainable fabric, comfortable for the baby, soft and flexible.

Non-elastic.

Pre-shaped for easy placing of the infant.

With strap for suspension from spring scale.

Seams with double stitching.

Trousers are able to hold infants up to 25 kg.

Dimensions of the trousers:

Height at centre: 27 cm

Waist: 29cm (58 cm circumference)

Side: 16 cm

Leg openings between 13 - 14cm (26 - 29cm circumference)

Crotch width between 11 - 15 cm

Suspension strap between 124 - 130 cm length, extending from leg opening along the side seams of trousers, fixed on both sides of the trousers.

---

Carry bag:

With shoulder strap, overlapping flap and closure, fits 5 weighing trousers.



**Laboratory Instrument India**