



Email : sales@laboratoryinstrumentindia.com

Product Name :
Weighing sling,infant,5kg

Product Code :
LBNY-0022-2200029



www.laboratoryinstrumentindia.com

Description :

One size fits all sling for spring scales, for newborn and infants aged 0-6 months.
Strong, durable and easy to clean fabric, comfortable for the baby, soft and flexible but non-elastic.
Seams with double stitching.
Square shape.
Lengths of the straps: 40cm, +/- 10%
Measurements of the sling: square shape, 80 x 80cm, +/- 5%
Can hold an infant up to 5 kg
With straps/handles for suspension from spring scale: 2 straps, each strap to be fitted with both ends on 2 parallel corners of the sling.

Material:

Material does not shrink more than 10% with the above indicated washing.
Machine washable up to 40°C with domestic type detergent.

Technical Specification :

WEIGHT AND VOLUME:

Volume: 0.0006m³

Weight: 0.1kg

Laboratory Instrument India

