



Email : sales@laboratoryinstrumentindia.com

Product Name :
Port. Photometer, Vitamin A in Veg. Oil

Product Code :
LBNY-0022-2200036



www.laboratoryinstrumentindia.com

Description :

Portable photometer to determine quantity of vitamin A content in refined edible oils. Measuring device equipped with light emitting diode (LED) technology and graphic LCD display.

Technical Specification :

Measuring device equipped with light emitting diode (LED) technology and graphic LCD display.

Light source: LED lamp

Usage: Portable

Staff qualification: 1 day training

Weight: Device: 0.45kg, case incl. accessories
and device:1.8kg

Delivered with accessories, assembled and packed in a portable case.

No need for user calibration.

Life span of approx. 30 000 hours of the LED - 10 million measurements.



Laboratory Instrument India

