



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Portable Fluorometer, Vitamin A in Foods

**Product Code :**

LBNY-0022-2200037



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Portable Fluorometer to determine quantity of Vitamin A content in foods and biological fluids. Measuring device equipped with light emitting diode (LED) technology and graphic LCD display.

**Technical Specification :**

Measuring device equipped with light emitting diode (LED) technology and graphic LCD display. Delivered with accessories. Assembled and packed in a portable case.

Light source: LED lamp

Use: Portable, point of care (POC)

Staff qualification: 1 day training

No need for user calibration.

Life span of approx. 30 000 hours of the LED - 10 million measurements.

Weight: Device: 0.45kg, case incl. accessories  
and device:1.9kg



**Laboratory Instrument India**

