



Email : sales@laboratoryinstrumentindia.com

Product Name :

Portable Fluorometer, Vitamin A in Foods

Product Code :

LBNY-0022-2200037



www.laboratoryinstrumentindia.com

Description :

Portable Fluorometer to determine quantity of Vitamin A content in foods and biological fluids. Measuring device equipped with light emitting diode (LED) technology and graphic LCD display.

Technical Specification :

Measuring device equipped with light emitting diode (LED) technology and graphic LCD display. Delivered with accessories. Assembled and packed in a portable case.

Light source: LED lamp

Use: Portable, point of care (POC)

Staff qualification: 1 day training

No need for user calibration.

Life span of approx. 30 000 hours of the LED - 10 million measurements.

Weight: Device: 0.45kg, case incl. accessories
and device:1.9kg



Laboratory Instrument India

