



Email : sales@laboratoryinstrumentindia.com

Product Name :
Portable Photometer, Iodine in Salt

Product Code :
LBNY-0022-220038



www.laboratoryinstrumentindia.com

Description :

Measuring device equipped with light emitting diode (LED) technology and graphic LCD display. Portable photometer to determine quantity of iodine content in salt. Assembled and packed in a portable case. Delivered with accessories.

Technical Specification :

Light source: LED lamp
Staff qualification: 1 day training
Use: Portable, point of care (POC)
No need for user calibration.
Life span of approx. 30 000 hours of the LED - 10 million measurements.
Weight: Device: 0.45kg, case incl. accessories and device: 1.9kg



Laboratory Instrument India