



Email : sales@laboratoryinstrumentindia.com

Product Name :
Portable Photometer, Iron in Foods

Product Code :
LBNY-0022-2200039



www.laboratoryinstrumentindia.com

Description :

Measuring device equipped with light emitting diode (LED) technology and graphic LCD display. Portable Fluorometer to determine quantity of iron in food. Delivered with accessories. Assembled and packed in a portable case.

Technical Specification :

Use: Portable, point of care (POC)
Light source: LED lamp
Staff qualification: 1 day training
Weight: Device: 0.45kg, case incl. accessories and device: 1.9kg



Laboratory Instrument India