



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Salt test kit, salt fort. w/IODIDE

**Product Code :**  
LBNY-0022-2200040



www.laboratoryinstrumentindia.com

**Description :**

Kit indicates presence or absence of iodine (if the salt is fortified with potassium iodide). Iodised salt test kit, for determining presence or absence of iodine in salt fortified with potassium iodide (KI).

**Technical Specification :**

Kit contains:

- a. Test solutions, two 10 ml ampoules/vials with white caps.
- b. A device (safety pin) for opening the ampoules or vials.
- c. Salt holder (standardized plastic measuring sample).
- d. Recheck solution for identifying false negatives (for alkaline samples), one red 10 ml ampoule / vial.
- e. Colour reference chart.
- f. Instruction leaflet(s) in English, French, Spanish and Arabic. Instruction leaflet in other languages can be provided upon request.



**Laboratory Instrument India**