



Email : sales@laboratoryinstrumentindia.com

Product Name :
Work Ware

Product Code :
LBNY-0004-820053



Description :

Work Ware

Technical Specification :

Work ware

Size : 16", Material: Leather

Leather welding gloves are guaranteed to with stand extreme temperafures up to 932 of (500) constructed with 100% soft sweat

Absorbent cotton lining - 100% soft sweat absorbent cotton lining for added softness and heat resistance

Ouzong Heat Resistant Leather Weltling Sleeves with Elastic Cuff Welding Apparel for Welding

Glass Transport and Processing 18" Golden

Brown Coal Forge Glove

Qeelink Leather Welding Apron - Heat & Flame Resistant Heavy Duly Work Apron with 6 Pockets, 42" Extra Large & Cross Back Extra Long Strap

3m Peltor Optime 98 Cap-mouut Earmuffs H9p3e0



Laboratory Instrument India

