



Email : sales@laboratoryinstrumentindia.com

Product Name :

Salt test kit, salt fort. w/IODATE

Product Code :

LBNY-0022-2200041



www.laboratoryinstrumentindia.com

Description :

Kit indicates presence or absence of iodine (if the salt is fortified with potassium iodate). Iodised salt test kit, for determining presence or absence of iodine in salt fortified with potassium iodate (KIO₃).

Technical Specification :

Kit contains:

- a. Test solutions, two 10 ml ampoules/vials with white caps.
- b. A device (safety pin) for opening the ampoules or vials.
- c. Instruction leaflet(s) in English, French, Spanish and Arabic. Instruction leaflet in other languages can be provided upon request.
- d. Salt holder (standardized plastic measuring sample).
- e. Color reference chart.
- f. Recheck solution for identifying false negatives (for alkaline samples), one red 10 ml ampoule / vial.



Laboratory Instrument India