



Email : sales@laboratoryinstrumentindia.com

Product Name :

Scale,survey,mother/child,bat,150kg

Product Code :

LBNY-0022-2200043



www.laboratoryinstrumentindia.com

Description :

Powered with customer replaceable batteries. The scale allows a child's weight to be measured while being held by an adult. Digital weighing scale to weigh adults and children, weighing capacity up to 150kg. with 100g. graduation.

Technical Specification :

Measuring range 0.1 to 150 kg.

Graduation: 100 g.

Accuracy: ± 100 g over the entire range of measurement.

Taring: Taring function is automatically initiated, without pressing any operational button. Taring function is initiated once a parent holds a child after being weighed first (mother child mode).

"Mother-child mode" is the default mode, alternative mother only mode can be activated by pressing an operational button.

Portability: Easy to carry by a health worker in the field. Maximum weight 4.0 kg (including batteries); maximum dimensions: length: 360 mm, width: 400 mm, height: 70 mm.

Power saving feature automatically turns off the scale after 4-6 minutes of non-use.

Main on/off button has a feature to prevent accidental switching on or off and is easily accessible.

The on/off button does not serve any other purpose and is encased to be protected from rain, dust and humidity

Light-coloured foot mat to avoid heating if exposed to direct sun light.

Footprint markings on the mat to indicate correct standing position (optional).

Display

Back-lit LCD or other feature for low light reading (eg 1600 lumens).

Fixed weight display occurs within 5 seconds and is maintained for a minimum of 20 seconds after the load is removed or until the next weighing.



Laboratory Instrument India