



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Scale, infant, tray, 20kgx5g, batteries

**Product Code :**

LBNY-0022-2200044



www.laboratoryinstrumentindia.com

**Description :**

Functions: automatic switch-off, TARE, HOLD, auto-calibration with each switch-on.

Large LCD display

Reading time max 5 seconds

Splash proof and shock resistant light-weight body

Smooth surface/finishing for easy cleaning/disinfection.

Weighing range 0 - 10kg

Minimum graduation: 5g.

Accuracy:  $\pm 5g$ .

Display in kilograms

All vital parts made of rust proof, child-friendly materials; display cover made of durable, transparent plastic.

Additional power supply: customer replaceable batteries.

Colour: white or similar: off-white, cream, light grey etc.

Power requirements: 220 V/50Hz

**Technical Specification :****ENVIRONMENTAL CONDITIONS:**

Storage conditions: Temperature:  $-10^{\circ}\text{C}$  -  $+65^{\circ}\text{C}$ , Humidity: 0-95%

Operating conditions: Temperature:  $10^{\circ}\text{C}$  -  $+40^{\circ}\text{C}$ , Humidity: 0-95%

**WEIGHT AND VOLUME:**

Volume:  $0.019\text{m}^3$

Weight: 2.60kg



**Laboratory Instrument India**