



Email : sales@laboratoryinstrumentindia.com

Product Name :

Scale, infant, tray, 20kgx5g, batteries

Product Code :

LBNY-0022-2200044



www.laboratoryinstrumentindia.com

Description :

Functions: automatic switch-off, TARE, HOLD, auto-calibration with each switch-on.

Large LCD display

Reading time max 5 seconds

Splash proof and shock resistant light-weight body

Smooth surface/finishing for easy cleaning/disinfection.

Weighing range 0 - 10kg

Minimum graduation: 5g.

Accuracy: $\pm 5g$.

Display in kilograms

All vital parts made of rust proof, child-friendly materials; display cover made of durable, transparent plastic.

Additional power supply: customer replaceable batteries.

Colour: white or similar: off-white, cream, light grey etc.

Power requirements: 220 V/50Hz

Technical Specification :

ENVIRONMENTAL CONDITIONS:

Storage conditions: Temperature: -10°C - $+65^{\circ}\text{C}$, Humidity: 0-95%

Operating conditions: Temperature: 10°C - $+40^{\circ}\text{C}$, Humidity: 0-95%

WEIGHT AND VOLUME:

Volume: 0.019m^3

Weight: 2.60kg



Laboratory Instrument India