



Email : sales@laboratoryinstrumentindia.com

Product Name :

MUAC,Adult,w/o color code,French

Product Code :

LBNY-0022-2200046



www.laboratoryinstrumentindia.com

Description :

Pack of 50 tapes with pictorial and text instructions for use printed on the reverse side, reusable. Adults' Mid Upper Arm Circumference (MUAC) measuring tape with a range of up to 50 cm.

Intended use:

To measure the mid upper arm circumference for rapid screening of nutritional status of pregnant/lactating women.

Technical Specification :

Tape for measuring adult's Mid Upper Arm Circumference (MUAC).

With three incisions to facilitate tape insertion and fixation, including one wider result window.

Print: Permanent, resistant to solvents; easily readable in low light working situations.

Double-sided, featuring user instructions in French ('tape method') on the reverse side.

Measuring range up to 50 cm.

Graduated with 1 mm precision.

Accuracy: ± 1 mm of the maximum measurement (50 cm).

Working/storage temperature range: 0 - 50 °C.

Material:

Materials used are free from phthalates, heavy metals and other dangerous matters as per the latest Candidate List of substances of very high concern for Authorisation published.

Non-tear stretch-resistant polypropylene.

Minimum thickness 0.3mm.



Laboratory Instrument India