



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
MUAC,Adult,w/o color code,Spanish

**Product Code :**  
LBNY-0022-220047



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

**Intended use:**

To measure the mid upper arm circumference for rapid screening of nutritional status of pregnant/lactating women.

Adults' Mid Upper Arm Circumference (MUAC) measuring tape with a range of up to 50 cm. Pack of 50 tapes with pictorial and text instructions for use printed on the reverse side, reusable.

**Technical Specification :**

Double-sided, featuring user instructions in Spanish('tape method') on the reverse side.

Print: Permanent, resistant to solvents; easily readable in low light working situations.

Tape for measuring adult's Mid Upper Arm Circumference (MUAC).

With three incisions to facilitate tape insertion and fixation, including one wider result window.

Working/storage temperature range: 0 - 50 °C.

Measuring range up to 50 cm.

Graduated with 1 mm precision.

Accuracy:  $\pm 1$  mm of the maximum measurement (50 cm).

**Material:**

Non-tear stretch-resistant polypropylene.

Minimum thickness 0.3mm.



**Laboratory Instrument India**