



Email : sales@laboratoryinstrumentindia.com

Product Name :
MUAC,Child 11.5 Red,PAC-50,Arabic

Product Code :
LBNY-0022-2200048



www.laboratoryinstrumentindia.com

Description :

To measure the mid upper arm circumference for rapid screening of nutritional status of children. Target population: Malnourished children or those at risk, aged 6-59 months.

Technical Specification :

Tape for measuring child's Mid Upper Arm Circumference (MUAC).
Double-sided, featuring user instructions in Arabic ('tape method') on the reverse side.
Print: Permanent, resistant to solvents; easily readable in low light working situations.
Working/storage temperature range: 0 - 50 °C.
Measuring range up to 26.5 cm.
Accuracy: ± 1 mm of the maximum measurement (26.5 cm).
Graduated with 1 mm precision with thicker line at 21.0 cm.
With three incisions to facilitate tape insertion and fixation, including one wider result window.

Material:

Minimum thickness 0.3 mm.
Non-tear stretch-resistant polypropylene.

Laboratory Instrument India

