



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
MUAC,Child 11.5 Red,PAC-50,Spanish

**Product Code :**  
LBNY-0022-2200050



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

To measure the mid upper arm circumference for rapid screening of nutritional status of children. Target population: Malnourished children or those at risk, aged 6-59 months.

**Technical Specification :**

Tape for measuring child's Mid Upper Arm Circumference (MUAC).  
Double-sided, featuring user instructions in Spanish ('tape method') on the reverse side.  
Print: Permanent, resistant to solvents; easily readable in low light working situations.  
With three incisions to facilitate tape insertion and fixation, including one wider result window.  
Working/storage temperature range: 0 - 50 °C.  
Measuring range up to 26.5 cm.  
Graduated with 1 mm precision with thicker line at 21.0 cm.  
Accuracy:  $\pm 1$  mm of the maximum measurement (26.5 cm).

**Material:**

Non-tear stretch-resistant polypropylene.  
Minimum thickness 0.3 mm.

---

# Laboratory Instrument India

