



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sachet,health card,160x220mm/BOX-500

Product Code :
LBNY-0023-220002



Description :

Material: Polyethylene

Sachet, transparent, open along 1 short side with a ziplock, securely sealed along the 2 long sides and one short side. The ziplock must fit 210 x 148mm(A5) size health cards.

Dimensions: length of 220mm, width of 160mm and thickness of 0.08mm (80 microns).

Technical Specification :

Plastic bag for health card, 160x220mm, Box of 500.



Laboratory Instrument India