



Email : sales@laboratoryinstrumentindia.com

Product Name :
Users guide School in a Box Kit, English

Product Code :
LBNY-0023-220006



Description :

Users guide School in a Box Kit, English Module 4

Technical Specification :

Paper: Inside, FSC certified 115 Matt Coated. Cover, FSC certified 250 GSM Matt Coated
Printing: Inside 4 colors CMYK recto/verso, Cover 4 colors recto/verso plus matt lamination recto only
Digital files will be submitted in high resolution PDF format, source files will also be included.
Size: Closed A4 (Portrait) Open A3
Number of pages: Inside 112, cover4
Binding: Spiral Binding (Wire O Wire)



Laboratory Instrument India