




Email : sales@laboratoryinstrumentindia.com

Product Name :
Users guide School in a Box Kit, French

Product Code :
LBNY-0023-220007



Description :

Users guide School in a Box Kit, French Module 4.

Technical Specification :

Printing: Inside 4 colors CMYK recto/verso, Cover 4 colors recto/verso plus matt lamination recto only
Digital files will be submitted in high resolution PDF format, source files will also be included.

Number of pages: Inside 122, cover4

Paper: Inside, FSC certified 115 Matt Coated. Cover, FSC certified 250

GSM Matt Coated

Size: Closed A4 (Portrait) Open A3

Binding: Spiral Binding (Wire O Wire).



Laboratory Instrument India