




Email : sales@laboratoryinstrumentindia.com

Product Name :
Trainer of Trainers Guide, Spanish

Product Code :
LBNY-0023-2200015



Description :

Size: Closed A4 (Portrait) Open A3
Number of pages: Inside 140, cover 4.
Binding: Spiral binding (Wire O Wire)
Printing: Inside 4 colors CMYK recto/verso, Cover 4 colors recto/verso plus matt lamination recto only
Paper: Inside, FSC certified, 115GSM, Matt Coated. Cover, FSC certified 250 GSM Matt Coated.

Technical Specification :

Trainer of Trainers Guide, Spanish



Laboratory Instrument India