




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Trainer of Trainers Guide, Spanish

**Product Code :**  
LBNY-0023-2200015



**Description :**

Size: Closed A4 (Portrait) Open A3

Number of pages: Inside 140, cover 4.

Binding: Spiral binding (Wire O Wire)

Printing: Inside 4 colors CMYK recto/verso, Cover 4 colors recto/verso plus matt lamination recto only

Paper: Inside, FSC certified, 115GSM, Matt Coated. Cover, FSC certified 250 GSM Matt Coated.

**Technical Specification :**

Trainer of Trainers Guide, Spanish



**Laboratory Instrument India**